

General Guide to the Linkam THMS600 Stages

APS Detector/Equipment Pool

Table of Contents

- Quick Guide
- Safety and Handling
- THMS600 Connections and Parts
- LN2 Sample Cooling
- Water Cooling for Stage Body
- Software Controls



Reference: much of the information in this guide is gleaned from User Manuals provided by Linkam with purchase of equipment. If more information is needed, please contact Detector Pool staff at 2-9493 or dp@aps.anl.gov.

The Linkam THMS600 stage can be used to heat or cool a sample from -196°C to 600°C at speeds of $1\text{-}100^{\circ}\text{C}$ per minute. The Detector Pool has one older stage as well as a newer stage with an updated controller design.

Safety & Handling

- Linkam stages have **EXTREMELY** fragile platinum leads; please do not touch or move.
- When heating above 300°C, stages require water cooling. Contact DP staff (2-9490).
- Regarding gas purging:
 - Do not use hydrogen or helium (thermal conductivity)
 - Dilute any reactive gases with 95% inert gas
 - Use only mixtures containing less than 20% oxygen
 - Flow less than 60 CC/min
- If using LN2 sample cooling, please follow all APS procedures for safe handling of LN2.



- The black capillary tube on the LN2 dewar lid is fragile; please handle carefully.
- Disconnect the LN2 dewar from the stage before heating above 300°C.

Hardware Operation


- Make all connections



TS1500



THMS600



Water cooling



DSC with cooling






- If cooling sample with LN2:
 - Fill the dewar approximately 2/3 full.
 - The LNP95 must be switched on before the T95/CI95 system controller.
 - The stage chamber needs to be purged of air before starting a cooling experiment.

Quick Guide - 2

Software Operation

- Turn on the computer, and log in
 - Username = dpuser
 - Contact DP staff for the password
 - Or, use the sector's LDAP account
- The start-up screen (lower left) will appear
 - Select "Linkam" from the dropdown menu
 - Click "Start" to start the IOC and medm
- Use the EPICS control screen (lower right) to set rate ($^{\circ}\text{C}/\text{min}$) and desired temperature.
 - Pump speed is generally set to "Automatic."
 - To adjust the y-axis on the strip chart display, right-click, then select "PV Limits."

Sample Loading




Aluminum and sapphire crucibles

Tzero Press, often used with DSC


Safety and Handling

- The heating element may be fragile. Exercise care when loading and unloading samples.
- Use as little sample as possible to reduce thermal load and avoid damaging the heating element.
- Heating over 300°C
 - Stages require **water cooling** when heated above 300°C. Contact DP staff (2-9490 or dp@aps.anl.gov) to obtain water circulator.
 - Disconnect the LN2 dewar from the stage before heating above 300°C. Leaving the tubing connected can damage the tubing and/or heater.
- Regarding gas purging:
 - Do not use hydrogen or helium (due to high thermal conductivity)
 - Do not use a gas mixture containing more than 20% oxygen
 - Dilute any reactive gases with 95% inert gas such as nitrogen or argon
 - Use a gas flow of less than 60 CC/min to avoid overloading the heating element
- Failure to observe these guidelines may result in heating element failure.
 - If you plan to use the stage in a manner not recommended above, please purchase your own heating element.
- If using LN2 sample cooling, please follow all APS procedures for safe handling of LN2.

Connections for THMS600 (old) with & w/o cooling



Crossover/null modem cable connects to computer serial port.




Tube connects to front of LNP. Make sure tube is secure.


PC cable connects CI94 to LNP

LN2 dewar connects to stage; securely connect inner tube first.

Connections for THMS600 (new) with & w/o cooling



Crossover/null modem cable connects to computer serial port.



Tube connects to front of LNP. Make sure tube is secure.


I²C cable connects CI94 to LNP

LN2 dewar connects to stage; securely connect inner tube first.

THMS600 Parts


Lid Window Assembly

To replace the windows in the Stage Lid (1) use the Window Tool (2) and align the two wide spacing pins to the Tube Clip Holder holes and unscrew the Lid Insert (3). The Stage Lid and Lid Insert should be turned upside down as shown in the diagram opposite and reassembled in the order indicated. The Lid Insert should be screwed down until the cover slips are held firmly, then turn the assembly over and screw down the Lid Insert until it is felt to come to a stop.



Bottom Window Assembly



Use two narrow spacing pins of Window Tool (2) to align it to the two holes of Window Locking Ring (4) and unscrew. Reassemble the bottom window as shown in the opposite diagram



- 15. Heating element wire
- 16. Stainless steel cooling tube
- 17. 22mm diameter pure silver heating block
- 18. Platinum temperature sensor
- 19. Sample holder ramp
- 20. Aperture hole
- 21. Vertical Sample Holder

Stage Assembly

1. Lemo connector for stage lead
2. Heating element carrier assembly
3. Stage body
4. Stage body water connector
5. Gas purge valve
6. Y-Sample manipulator
7. Stage door
8. Door locking thumbscrew
9. X-Sample manipulator
10. Liquid nitrogen cooling connector
11. Bypass stage body water cooling connector
12. Sample chamber




Sample Cooling with LNP95

- If using LN2 cooling, please follow all APS procedures for safe handling of LN2.
- The black capillary tube on the dewar lid is fragile. Please handle carefully.
- Fill the dewar approximately 2/3 full.
- The LNP95 must be switched on before the T95/CI95 system controller.
- The stage chamber needs to be purged of air before starting a cooling experiment.
 - See purging procedure on next page
- Disconnect the LN2 dewar from the stage before heating above 300°C. Leaving the tubing connected can damage the tube and/or heater.



Purging Procedure

- For full details, please see the THMS600 User Guide from Linkam, available on the Detector Pool website.
- To purge using recycled nitrogen gas produced from the Dewar:
 - Make all connections, make sure that the stage lid and gas inserts are secure, turn on the LNP95, then turn on the T95/CI95.
 - Set the LNP95 to “Manual” mode.
 - Set the T95/CI95 to hold at 40°C.
 - Set the LNP95 to a speed of 100.
 - Block the hole in the white plastic pump connector, and pinch the narrow window tube to block it.
 - Block the gas outlet for a few seconds to allow pressure to build, then release. Repeat for a few minutes.
 - When finished, unblock the pump connector, window, tube, and gas outlets.
 - Change the LNP95 to “Automatic” mode.




Pinch window tube
with left hand

Block hole in
Connector with finger
pump connector


Purging
tube

Block and release outlet
valve with finger

Water Cooling Connections




- When heating above 300°C, stages require water cooling to keep the stage body and lid cool. Contact DP staff (2-9490; dp@aps.anl.gov) for a water circulator.
- There are several options for cooling:
 - Linkam ECP Water Circulating Pump
 - It may be necessary to prime the ECP if the water connectors have been removed.
 - DP circulating pump with custom housing
 - To ensure proper circulation, be careful to avoid introducing air into the system. While connecting and disconnecting, take care to prevent water from leaking out.
 - Beamlines may use their own small, recirculating chiller. (less than 20CC/min flow rate)



EPICS Software Controls

- Turn on the computer, and log in.
 - Username = dpuser
 - Contact DP staff for the password
 - Or, use the sector's LDAP account (see beamline personnel)
- The start-up screen (shown right) will appear.
 - Select "Linkam" from the dropdown menu
 - Click "Start" to start the IOC and medm
- The figure below shows the two resulting terminal windows and the Linkam control screen.



- The following parameters can be set (left side):
 - Pump speed for LN2 cooling (optional—this can also be controlled automatically based on the selected cooling rate). Upper limit = 30
 - Rate of temperature change (°C/min)
 - Desired temperature
- Note the control buttons (right side):
 - Scan Disable: tells VME (EPICS sscan) to stop communicating with the controller (e.g. software is loaded but controller is not in use). "Enable" restarts communication.
 - Toggle LNP between automatic & manual
 - Start, hold, or stop heating/cooling

EPICS Software Controls

ci94_new.adl

Linkam CI94

Error Status

- Cooling Rate too fast
- Open Circuit
- Power Surge
- No Exit 300
- Both Stages
- Link Error
- NC
- Default

Scan Disable
Enabled

Pump Control
Automatic

	Set Value	Current Value
Pump Speed	0	0
Rate	0.00	
Temperature 1	0.0	28.00
Temperature 2		28.00
DSC value		0
Status	Stopped	

Start
Hold
Stop
Clear Buffer
End After Profile
End On Stop

Temperature 1

Temp (c)

Print
Close
PV Info
PV Limits
MEDM Main Window
Display List
Toggle Hidden Button Markers
Refresh
Retry Connections

Strip Chart Data

Channel	Color	Low Source	Low Limit	Channel Default	High Limit
xp_linkam_xrd45:ci94:temp	Blue	User Specified	0.00	User Specified	1200.00
xp_linkam_xrd45:ci94:dsc	Purple	Channel	0.00	Channel	1.00
...		Channel	0.00	Channel	1.00
...		Channel	0.00	Channel	1.00
...		Channel	0.00	Channel	1.00
...		Channel	0.00	Channel	1.00
...		Channel	0.00	Channel	1.00
...		Channel	0.00	Channel	1.00
...		Channel	0.00	Channel	1.00
...		Channel	0.00	Channel	1.00

Period: 10 Units: minute

Apply Cancel Help



To adjust the scale on the Y-axis: right-click, then select "PV Limits."

Select "User Specified," and set the high and low limits. Click "Apply."

Thanks to John Hammonds for the EPICS interface!

Appendix: Linksys Software Controls

- The Linkam equipment can be controlled with either EPICS or vendor software called Linksys, but not both simultaneously.
 - The Detector Pool no longer supports the Linksys software, but some beamlines may have it installed on their computers.
 - To use the Linkam software, make sure that the EPICS IOC is not running.
- From the desktop, click Linksys32 icon
- Select “File: Connect”




Start and stop temperature ramps

Set the rate, max/min temperature, and hold time here

Data capture controls

Appendix: Linksys Software Controls

Select "View Temperature Profile" to setup multiple-step temperature ramps



A previously-set temperature profile may affect your ability to control the stage. If control difficulties occur, check the temperature profile settings.

Profile				
Profile - Cycle mode off				
Ramp	Rate	Limit	Time	Delay
1	1	20.0	30	-
2	0	0.0	0	-
3	0	0.0	0	-
4	0	0.0	0	-
5	0	0.0	0	-
6	0	0.0	0	-
7	0	0.0	0	-
8	0	0.0	0	-
9	0	0.0	0	-
10	0	0.0	0	-
11	0	0.0	0	-

If you are using the DSC, you may be required to use a Temperature Profile!!!